



ACCEL® III


Use Instructions



L'électronique au service de l'automobile

Table of contents

PRESENTATION.....	4
USE	5
TECHNICAL CHARACTERISTICS	7



ACCEL® III Presentation

ACCEL® III consists of two identical acceleration commands positioned on both sides of the steering wheel. Pushing down of these levers allows for a very precise acceleration.



Placed under the steering wheel, at the level of the column casing, it permits to keep the airbag.

Two possibilities for driving :

- conventional, with accelerator pedal
- manual, with *ACCEL® III*

An electronic push button with a warning light activates the manual mode.



ACCEL® III can be installed on vehicles equipped with automatic gear and electrical accelerator pedal.

A *Brake on steering wheel* should be installed as a complement to this system.

ACCEL® III Use

1- At the moment of starting up, the conventional mode of driving is active.

2- Selection of manual mode :

- A 30 second pressure on the push button after the starting up of the engine activates the manual mode, warning light turned on to confirm this selection.




It is impossible to interrupt manual mode without switching off the engine.

- **If manual mode is not selected after 30 seconds, it is necessary to restart the engine and repeat the selection procedure.**

Security - Applying the brakes stops acceleration when manual mode is activated.



ACCEL® III makes possible a very precise acceleration. Starting up and manoeuvring are easier and full power of the vehicle is maintained.



Acceleration is progressive at the beginning of the levers travel.

Technical characteristics of **ACCEL® III**

Working temperature :

- from – 40° to + 85°C

Power Supply :

- 12 V (from 10 to 16 V)

Complies with regulation 95/54/EC





SOJADIS EQUIPEMENT
Place du Maréchal Leclerc – 49 510 JALLAIS – France
S.A.R.L. au capital de 7500€ RCS ANGERS 2007 B 00057
Siret 493 658 785 00016 - Code APE 4520 B - TVA intra : FR58 493 658 785



Contact : +33 (0)2 41 64 05 58
sojadis@sojadis.com